



Product list 2020



PharmaMedia Dr. Müller GmbH (PMM) was founded in January 2013 by Dr. Rolf Müller, who has also founded heipha Dr. Müller GmbH more than 40 years ago.

The focus of **PMM** is the production of high-quality culture media as well as supplying state of the art service to customers of the pharmaceutical industry. The production facility in Leimen, Germany, is especially designed to assure the highest level of security.

We have succeeded in manufacturing innovative products which offer a high level of quality and product security. Our product portfolio extends continuously.

We will be pleased to be at your disposal for any question concerning culture media and we are looking forward to receiving your requests.

PMM GmbH

Gustav-Throm-Str. 1

69181 Leimen

Phone: +49(0) 6224 92807-0

Orders: order@pmm-leimen.de

E-Mail: info@pmm-leimen.de

www.pmm-leimen.de

Environmental Monitoring (Cleanrooms class A and B as well as in isolators)

CSG- Click& Safe gamma

- triple bagged, gamma-irradiated
- lockable petridishes
- neutralizer: LTHT (L = lecithin/T = tween/H = histidine/T = Na-thiosulphate/
G = Glycine)
- long shelf-life
- room temperature storage: 15-25°C
- temperature controlled delivery on request
- long incubations possible
- H₂O² – impermeable foil
- unique identification assigned by each plate with label and barcode

Contact plates			
article number	description	unit	
100.0100	TSA+ LTHT	100	
120.0100	SDA + LTHT	100	
101.0100	TSA U+ To neutralise disinfectants (also for quaternary ammonium compound and biguanide)	100	
114.0100	TSA +LTG +β-Lactamase 2G + (incl. Penase) β-LactamasePlus: For the inactivation of penicillines as well as first to fifth generation of cephalosporins	100	
116.0100	TSA +LT + β-Lakt I/II	100	

Settle plates			
article number	description	unit	
200.0060	TSA + LTHT	60	
203.0060	TSA U+	60	On request
210.0060	TSA	60	
220.0060	SDA + LTHT	60	
214.0060	TSA +LTG + β -Lactamase 2G + (incl. Penase) β-LactamasePlus: For the inactivation of penicillines as well as first to fifth generation of cephalosporins	60	
216.0060	TSA +LT + β -Lakt I/II	60	
230.0060	R2A (gamma-irradiated, single-bagged)	60	On request

Transport bags			
article number	description	unit	
880.0120	Gamma-Bag	120	

Environmental Monitoring CS – Click & Safe

- single bagged
- lockable petridishes
- neutralizer: LTH (L = lecithin/T = tween/H = histidine)
- long shelf-life
- room temperature storage: 15-25°C
- temperature controlled delivery on request
- long incubations possible
- unique identification assigned by each plate with label and barcode

Contact plates			
article number	description	unit	
300.0120	TSA + LTH	120	
120.0100	SDA + LTHT (CSG)	120	

Settle plates			
article number	description	unit	
400.0060	TSA + LTH	60	
410.0060	TSA	60	
425.0060	SDA + Chloramphenicol	60	
426.0060	SDA	60	

Selective culture media CS - Click&Safe

- single bagged
- lockable petridishes
- neutralizer: LTH (L = lecithin/T = tween/H = histidine)
- long shelf-life
- room temperature storage: 15-25°C
- long incubations possible
- if request, bulk buying possible
- unique identification assigned by each plate with label and barcode

90 mm plates			
article number	description	unit	
400.0060	TSA + LTH	60	
410.0060	TSA	60	
425.0060	SDA + Chloramphenicol	60	
426.0060	SDA	60	
450.0060	Cetrimide	60	
455.0060	Mannitol Salt	60	
460.0060	McConkey	60	
465.0060	VRBD	60	
470.0060	XLD	60	
490.0060	R2A	60	
230.0060	R2A (gamma-irradiated, single-bagged)	60	On request

Culture media in bottles

- shrink-wrapped – packaging concept
- innovative bottle/cap-concept
- indicator labels for a safe autoclaving detection
- unique identification assigned by each plate with label and barcode

Volumes

- 100 ml infusion bottle (B=screw cap/ S=punctured plug)
- 200 ml screw cap
- 500 ml infusion bottle (B=screw cap/ S=punctured plug)
- 1 l screw cap

Screw-cap bottles			
article number	description	volumes	bottle
500.B100	TSB – Tryptic Soy Broth	100 ml	1
505.D000	TSB	1000 ml	1
542.D000	Fluid A	1000 ml	1
571.B090	NBP	90 ml	1
578.D000	NBP	1000 ml	1
580.B100	Sabouraud-Dextrose-Broth	100 ml	1
581.B100	McConkey Bouillon	100 ml	1
600.B200	TSA	200 ml	1
601.B200	SDA	200 ml	1
602.B200	SDA + Chloramphenicol	200 ml	1
603.B200	R2A	200 ml	1

Infusion bottles, Septum, Coloured Flip Caps			
article number		volumes	bottle
520.S100	TSB – Tryptic Soy Broth	100 ml	1
530.S100	FTM – Thioglycollate Clear	100 ml	1
540.S100	Fluid A	100 ml	1
550.S100	Fluid D	100 ml	1

We are pleased to offer you special productions!!

Information and order

Email: order@pmm-leimen.de

Tel.: +49 (0)6224-928070

Fax: +49 (0)6224-9280769

Terms and conditions: www.pmm-leimen.de/downloads/

➔ AGB