



Environmental Monitoring (Cleanrooms class A and B as well as in isolators) - CSG- Click& Safe gamma

- triple bagged, gamma-irradiated
- lockable petri dishes
- neutralizers: LTHT (L = lecithin / T = Tween 80 / H = histidine / T = Na-thiosulfate / G = Glycine)
- long shelf-life
- room temperature storage: 15-25°C
- temperature-controlled delivery on request
- long incubations possible
- H²O² – impermeable foil
- unique identification assigned by each plate with label and barcode

Contact plates			
Article number	Description	Units per pack	
100.0100	TSA+ LTHT	100	
120.0100	SDA + LTHT	100	
101.0100	TSA U+	100	
	To neutralize disinfectants (also for quaternary ammonium compounds and biguanides)		
119.0100	TSA +LTG +β-Lactamase + (incl. Penase)	100	
	β-LactamasePlus: For the inactivation of penicillins as well as first to fifth generation of cephalosporins		
Settle plates			
200.0060	TSA + LTHT	60	
210.0060	TSA	60	
220.0060	SDA + LTHT	60	
219.0060	TSA +LTG +β-Lactamase+ (incl. Penase)	60	
	β-LactamasePlus: For the inactivation of penicillins as well as first to fifth generation of cephalosporins		
230.0060	R2A (gamma-irradiated, single-bagged)	60	
Transport bags			
880.0120	Gamma-Bag	120	

Environmental Monitoring CS – Click & Safe

- single bagged
- lockable petri dishes
- neutralizer: LTH (L = lecithin/T = tween/H = histidine)
- long shelf-life
- room temperature storage: 15-25°C
- temperature controlled delivery on request
- long incubations possible
- unique identification assigned by each plate with label and barcode

Contact plates			
Article number	Description	Units per pack	
300.0120	TSA + LTH	120	

Settle plates			
Article number	Description	Units per pack	
400.0060	TSA + LTH	60	
410.0060	TSA	60	
426.0060	SDA	60	

Selective culture media CS - Click&Safe

- single bagged
- lockable petridishes
- neutralizer: LTH (L = lecithin/T = tween/H = histidine)
- long shelf-life
- room temperature storage: 15-25°C
- long incubations possible
- if request, bulk buying possible
- unique identification assigned by each plate with label and barcode

90 mm plates			
Article number	Description	Units per pack	
400.0060	TSA + LTH	60	
410.0060	TSA	60	
425.0060	SDA + Chloramphenicol	60	
426.0060	SDA	60	
450.0060	Cetrimide	60	
455.0060	Mannitol Salt	60	
460.0060	McConkey	60	
465.0060	VRBD	60	
470.0060	XLD	60	
490.0060	R2A	60	
230.0060	R2A (gamma-irradiated, single-bagged)	60	

Culture media in bottles

- shrink-wrapped – packaging concept
- innovative bottle/cap-concept
- indicator labels for a safe autoclaving detection
- unique identification assigned by each plate with label and barcode

Volumes

- 100 ml infusion bottle (B=screw cap/ S=punctured plug)
- 200 ml screw cap
- 500 ml infusion bottle (B=screw cap/ S=punctured plug)
- 1 l screw cap

Screw-cap bottles			
Article number	Description	Volume	Units per pack
500.B100	TSB – Tryptic Soy Broth	100 ml	12
505.D000	TSB	1000 ml	6
542.D000	Fluid A	1000 ml	6
552.D000	Fluid D	1000 ml	6
571.B090	NBP	90 ml	12
578.D000	NBP	1000 ml	6
580.B100	Sabouraud-Dextrose-Broth	100 ml	12
581.B100	McConkey Bouillon	100 ml	12
600.B200	TSA	200 ml	12
601.B200	SDA	200 ml	12
602.B200	SDA + Chloramphenicol	200 ml	12
603.B200	R2A	200 ml	12

Infusion bottles, Septum, Colored Flip Caps			
Article number	Description	Volume	Units per pack
520.S100	TSB – Tryptic Soy Broth	100 ml	20
530.S100	FTM – Thioglycollate Clear	100 ml	20
540.S100	Fluid A	100 ml	20
550.S100	Fluid D	100 ml	20

We are pleased to offer you special productions!!

Information and order

Email: order@pmm-leimen.de

Tel.: +49 (0)6224-928070

Fax: +49 (0)6224-9280769

Terms and conditions: www.pmm-leimen.de/download/

➔ AGB