



# Productlist



**PharmaMedia Dr. Müller GmbH (PMM)** was founded in January 2013 by Dr. Rolf Müller, who has also founded heipha Dr. Müller GmbH more than 40 years ago.

The focus of **PMM** is the production of high-quality culture media as well as supplying state of the art service to customers of the pharmaceutical industry. The production facility in Leimen, Germany, is especially designed to assure the highest level of security.

We have succeeded in manufacturing innovative products which offer a high level of quality and product security. Our product portfolio extends continuously.

We will be pleased to be at your disposal for any question concerning culture media and we are looking forward to receiving your requests.

**PMM GmbH**

Gustav-Throm-Str. 1

69181 Leimen

**Phone: +49(0) 6224 92807-0**

Orders: [order@pmm-leimen.de](mailto:order@pmm-leimen.de)

E-Mail: [info@pmm-leimen.de](mailto:info@pmm-leimen.de)

[www.pmm-leimen.de](http://www.pmm-leimen.de)

## Environmental Monitoring (Cleanrooms class A and B as well as in isolators)

### CSG- Click& Safe gamma

- triple bagged, gamma-irradiated
- lockable petri dishes
- neutralizer: LTHT (L = lecithin/T = tween/H = histidine/T = Na-thiosulfate/ = Glycine) G
- long shelf-life
- room temperature storage: 15-25°C
- temperature controlled delivery on request
- long incubations possible
- H<sup>2</sup>O<sup>2</sup> – impermeable foil
- unique identification assigned by each plate with label and barcode

Contact plates				
article number	description	unit	list price, net	
100.0100	TSA+ LTHT	100		
120.0100	SDA + LTHT	100		
101.0100	TSA U+ To neutralise disinfectants (also for quaternary ammonium compound and biguanides)	100		
119.0100	TSA +LTG +β-Lactamase + (incl. Penase) β-LactamasePlus: For the inactivation of penicillins as well as first to fifth generation of cephalosporins	100		

<b>Settle plates</b>				
<b>article number</b>	<b>description</b>	<b>unit</b>	<b>list price, net</b>	
200.0060	TSA + LTHT	60		
210.0060	TSA	60		
220.0060	SDA + LTHT	60		
219.0060	TSA +LTG + $\beta$ -Lactamase+ (incl. Penase) <b><math>\beta</math>-LactamasePlus:</b> For the inactivation of penicillins as well as first to fifth generation of cephalosporins	60		
230.0060	R2A (gamma-irradiated, single-bagged)	60		On request

<b>Transport bags</b>				
<b>article number</b>	<b>description</b>	<b>unit</b>	<b>list price, net</b>	
880.0120	Gamma-Bag	120		

## Environmental Monitoring CS – Click & Safe

- single bagged
- lockable petri dishes
- neutralizer: LTH (L = lecithin/T = tween/H = histidine)
- long shelf-life
- room temperature storage: 15-25°C
- temperature controlled delivery on request
- long incubations possible
- unique identification assigned by each plate with label and barcode

<b>Contact plates</b>			
<b>article number</b>	<b>description</b>	<b>unit</b>	<b>list price, net</b>
300.0120	TSA + LTH	120	
120.0100	SDA + LTHT (CSG)	120	

<b>Settle plates</b>			
<b>article number</b>	<b>description</b>	<b>unit</b>	<b>list price, net</b>
400.0060	TSA + LTH	60	
410.0060	TSA	60	
426.0060	SDA	60	

## Selective culture media CS - Click&Safe

- single bagged
- lockable petri dishes
- neutralizer: LTH (L = lecithin/T = tween/H = histidine)
- long shelf-life
- room temperature storage: 15-25°C
- long incubations possible
- if request, bulk buying possible
- unique identification assigned by each plate with label and barcode

<b>90 mm plates</b>			
<b>article number</b>	<b>description</b>	<b>unit</b>	<b>list price, net</b>
400.0060	TSA + LTH	60	
410.0060	TSA	60	
425.0060	SDA + Chloramphenicol	60	
426.0060	SDA	60	
450.0060	Cetrimide	60	
455.0060	Mannitol Salt	60	
460.0060	McConkey	60	
465.0060	VRBD	60	
470.0060	XLD	60	
490.0060	R2A	60	
230.0060	R2A (gamma-irradiated, single-bagged)	60	On request

## Culture media in bottles

- shrink-wrapped – packaging concept
- innovative bottle/cap-concept
- indicator labels for a safe autoclaving detection
- unique identification assigned by each plate with label and barcode

## Volumes

- 100 ml infusion bottle (B=screw cap/ S=punctured plug)
- 200 ml screw cap
- 500 ml infusion bottle (B=screw cap/ S=punctured plug)
- 1 l screw cap

<b>Screw-cap bottles</b>				
<b>article number</b>	<b>description</b>	<b>volumes</b>	<b>unit</b>	<b>list price, net</b>
500.B100	TSB – Tryptic Soy Broth	100 ml	12	
505.D000	TSB	1000 ml	6	
542.D000	Fluid A	1000 ml	6	
552.D000	Fluid D	1000 ml	6	
571.B090	NBP	90 ml	12	
578.D000	NBP	1000 ml	6	
580.B100	Sabouraud-Dextrose-Broth	100 ml	12	
581.B100	McConkey Bouillon	100 ml	12	
600.B200	TSA	200 ml	12	
601.B200	SDA	200 ml	12	
602.B200	SDA + Chloramphenicol	200 ml	12	On request
603.B200	R2A	200 ml	12	On request

<b>Infusion bottles, Septum, Coloured Flip Caps</b>			
<b>article number</b>		<b>volumes</b>	<b>unit list price, net</b>
520.S100	TSB – Tryptic Soy Broth	100 ml	20
530.S100	FTM – Thioglycollate Clear	100 ml	20
540.S100	Fluid A	100 ml	20
550.S100	Fluid D	100 ml	20 On request

**We are pleased to offer you special productions!!**

**Information and order**

Email: [order@pmm-leimen.de](mailto:order@pmm-leimen.de)

Tel.: +49 (0)6224-928070

Fax: +49 (0)6224-9280769

**Terms and conditions:** [www.pmm-leimen.de/angebot/](http://www.pmm-leimen.de/angebot/)

➔ AGB