

Technical Data Sheet													
Use in	<ul style="list-style-type: none"> Pharmaceutical Industry Cosmetic Industry Food Industry For industrial, laboratory & research applications only Basic medium according to EP 2.6.1 and USP <71> 												
Use for	<ul style="list-style-type: none"> Test for specified micro-organisms Preparation of test strains Diluent for sample preparation Identification and growth of aerobic micro-organisms 												
Typical composition per liter	<table> <tbody> <tr> <td>Casein peptone</td> <td>17g</td> <td>Glucose x H₂O</td> <td>2,5g</td> </tr> <tr> <td>Soy peptone</td> <td>3g</td> <td>K₂HPO₄</td> <td>2,5g</td> </tr> <tr> <td>NaCl</td> <td>5g</td> <td></td> <td></td> </tr> </tbody> </table> <p>This medium can be adjusted / or supplemented according to the performance criteria required.</p>	Casein peptone	17g	Glucose x H ₂ O	2,5g	Soy peptone	3g	K ₂ HPO ₄	2,5g	NaCl	5g		
Casein peptone	17g	Glucose x H ₂ O	2,5g										
Soy peptone	3g	K ₂ HPO ₄	2,5g										
NaCl	5g												
Filling volume	<ul style="list-style-type: none"> 1000 mL 												
Bottle format	<ul style="list-style-type: none"> 1100 mL laboratory glass screw cap bottle Type I glass (borosilicate glass) Bottle opening about 30mm Colour of cap: green GL45 screw cap with 3 integrated septa 												
Units per pack	<ul style="list-style-type: none"> 6 bottles on a plastic tray wrapped with shrink foil 												
Shelf life	<ul style="list-style-type: none"> 12 months from production date 												
Storage conditions	<ul style="list-style-type: none"> Recommended storage temperature: 2-25 °C Should be stored at temperatures as stable as possible Store protected from light exposure 												
Label	<ul style="list-style-type: none"> On the side Contain autoclave indicator 												
Label information	<ul style="list-style-type: none"> Product name: TSB 1000 mL Expiry date: YYYYMMDD → MMM in letters (e.g.: 2023Nov04) Lot-number Individual number Barcode 												
Barcode	<ul style="list-style-type: none"> 2-dimensional (data matrix), 20 digits: Digits 1-3: Art.-No. Digits 4-9: Lot-Number Digits 10-14: Individual-Number Digits 15-20: Date (YYMMDD) 												

Technical Data Sheet	
Delivery	<ul style="list-style-type: none"> • Temperature controlled delivery on request • For shipments of larger amounts plastic pallets in Euro-size can be used
Bottle information	<ul style="list-style-type: none"> • Label contains autoclaving indicator (brown → green) • Bottles are incubated for >48 hours at 25-35°C after production and then packed immediately • Bottles are not touched any more by hand after autoclaving
Place of production	PharmaMedia Dr. Müller GmbH Gustav-Throm-Str. 1, 69181 Leimen - Germany

Quality control, Certificates																																																										
Certificates	<p>Each lot of product can be obtained with a certificate of analysis (CoA):</p> <table border="1"> <thead> <tr> <th colspan="2">Physico-chemical test parameters:</th> </tr> </thead> <tbody> <tr> <td>Appearance</td> <td>Clear, yellowish</td> </tr> <tr> <td>pH value</td> <td>7,1 – 7,5</td> </tr> <tr> <td>Filling volume</td> <td>980 – 1040 mL</td> </tr> <tr> <td colspan="2"> </td> </tr> <tr> <th colspan="2">Growth Promotion test: 10-100 CFU</th> </tr> <tr> <td><i>S.aureus</i></td> <td>ATCC 6538</td> <td>30-35°C</td> <td>1 day</td> <td>Good growth</td> </tr> <tr> <td><i>E.coli</i></td> <td>ATCC 8739</td> <td>30-35°C</td> <td>1 day</td> <td>Good growth</td> </tr> <tr> <td><i>P.paraeruginosa</i></td> <td>ATCC 9027</td> <td>30-35°C</td> <td>1 day</td> <td>Good growth</td> </tr> <tr> <td><i>B.spizizenii</i></td> <td>ATCC 6633</td> <td>30-35°C</td> <td>1 day</td> <td>Good growth</td> </tr> <tr> <td><i>C.albicans</i></td> <td>ATCC 10231</td> <td>20-25°C</td> <td>3-5 d</td> <td>Good growth</td> </tr> <tr> <td><i>A.brasiliensis</i></td> <td>ATCC 16404</td> <td>20-25°C</td> <td>3-5 d</td> <td>Good growth</td> </tr> <tr> <td colspan="5"> </td> </tr> <tr> <th colspan="5">Sterility control</th> </tr> <tr> <td colspan="5">≥ 7 days at 30-35°C, no growth</td> </tr> </tbody> </table> <p>*In case of a direct inoculation and incubation in the bottle, please ensure that sufficient aeration of the bottle is warranted</p>	Physico-chemical test parameters:		Appearance	Clear, yellowish	pH value	7,1 – 7,5	Filling volume	980 – 1040 mL			Growth Promotion test: 10-100 CFU		<i>S.aureus</i>	ATCC 6538	30-35°C	1 day	Good growth	<i>E.coli</i>	ATCC 8739	30-35°C	1 day	Good growth	<i>P.paraeruginosa</i>	ATCC 9027	30-35°C	1 day	Good growth	<i>B.spizizenii</i>	ATCC 6633	30-35°C	1 day	Good growth	<i>C.albicans</i>	ATCC 10231	20-25°C	3-5 d	Good growth	<i>A.brasiliensis</i>	ATCC 16404	20-25°C	3-5 d	Good growth						Sterility control					≥ 7 days at 30-35°C, no growth				
Physico-chemical test parameters:																																																										
Appearance	Clear, yellowish																																																									
pH value	7,1 – 7,5																																																									
Filling volume	980 – 1040 mL																																																									
Growth Promotion test: 10-100 CFU																																																										
<i>S.aureus</i>	ATCC 6538	30-35°C	1 day	Good growth																																																						
<i>E.coli</i>	ATCC 8739	30-35°C	1 day	Good growth																																																						
<i>P.paraeruginosa</i>	ATCC 9027	30-35°C	1 day	Good growth																																																						
<i>B.spizizenii</i>	ATCC 6633	30-35°C	1 day	Good growth																																																						
<i>C.albicans</i>	ATCC 10231	20-25°C	3-5 d	Good growth																																																						
<i>A.brasiliensis</i>	ATCC 16404	20-25°C	3-5 d	Good growth																																																						
Sterility control																																																										
≥ 7 days at 30-35°C, no growth																																																										
Certificate of origin	<p>All media lots produced by PMM can be obtained with a Certificate of Origin (CoO). All animal derived raw materials are specified as follows:</p> <ul style="list-style-type: none"> • Raw material • Tissue • Animal source • Country of origin • Infectivity category (acc. to TSE guideline: EMA/410/01 rev. 3) 																																																									

Quality control, Certificates	
BSE policy	<ul style="list-style-type: none"> In compliance with the current note for guidance on minimizing the risk of transmitting animal spongiform encephalopathy via human or veterinary medicinal products, we check the CoO of raw material in respect to the specified animal source, the country of origin and the infectivity category. We neither store or process ruminant raw materials obtained from high infectivity tissues (IA) nor ruminant raw materials whose animal source originates from countries or regions with an undetermined risk (cat C/GBR IV).

Safety Data	
Toxic ingredients	<ul style="list-style-type: none"> None
Basic composition	<ul style="list-style-type: none"> See typical composition
Solvent content	<ul style="list-style-type: none"> None
Safety data sheet required	<ul style="list-style-type: none"> Not mandatorily required